

# Contents of the KariesScreenTest + P

The KariesScreenTest + P enables the determination of the bacteriological and functional saliva parameters and thus enables early detection of the bacterial shift from homeostasis to a dysbiotic state.

Bacterial detection for SM and LB



Buffer capacity test



Paraffin capsules to collect the saliva sample



CO<sup>2</sup> tablet



Measuring pipette 1 ml



Buffer solution



pH measuring strips to determine the buffer capacity




pH measuring strips for determining the saliva pH value



Calibration to measure the secretion rate

evaluation chart



Content:  
5 SM / LB  
5 buffer tests